

Endless Plus

“ Our most advanced lenses,
designed for a connected lifestyle.
Effortlessly clear and relaxed vision
whatever the day brings”



Spatial
Vision



Eye Focus
Profile



Ray
Tracing



Digital
Vision



Fitting
Configuration



For special, occasional use

**Easy to adapt.
Easy to wear.**

100% of wearers say they
feel **VERY SATISFIED**
with the lens performance



★ These principles underpin IOT's cutting-edge lens optimization technology,
IOT Digital Ray-Path 2 Technology.

iot

See the difference

Endless Plus

Endless Plus lenses are designed to alleviate the symptoms of eyestrain associated with the use of electronic devices. They incorporate a small increase in near power that can be customized for each patient.

- More relaxed vision requiring less accommodative effort
- Designed to significantly improve reading speed on digital devices
- Comfortable and precise focus at all distances
- Near elimination of peripheral blur
- Superior visual quality for viewing digital devices
- Impeccable visual quality

Ground breaking technology



Spatial
Vision

99.5%

of gaze directions are
now optimized

Use of accommodation for improved spatial vision and super dynamic vision.

“ Naturally good vision, even when you’re on the move”


Spatial Vision Technology optimizes the range within which an object can be moved closer to or farther away from the eye without a noticeable change in the focus of the image. This improves the eye’s ability to focus on objects both near and far, resulting in **exceptional sharpness perception and extended areas of focus.**

This technology ensures that wearers experience clear and sharp vision, even in dynamic environments where objects are constantly changing distances from the eye. It enables **smooth and rapid transitions across visual zones**, especially when looking down to focus on nearby objects.

**It enhances your spatial vision,
providing extended areas of clear focus.**



WITHOUT SPATIAL VISION TECHNOLOGY



WITH SPATIAL VISION TECHNOLOGY



Eye Focus
Profile

Eye Focus profiling for flawless vision.

“Flawless vision from any angle”

Eye Focus Profile Technology maximizes the clear viewing area and significantly reduces perceived blur in the peripheral areas of the lens, achieving an **absolute reduction in defocus and increasing the visual performance of the lenses.**



Ray
Tracing

Optimum vision across the entire lens.

“Unlock superior near vision precision perfect for your connected, on-the-go lifestyle”

Thanks to point-by-point optimization over the entire lens surface, Ray Tracing Technology provides **precise vision at every distance and in every direction of gaze.** Wearers enjoy a more natural visual experience.

Optimum vision across the entire lens



Digital
Vision

Comfortable focus when using digital devices.

“Embrace a connected, on-the-go lifestyle”

These advanced lenses feature a **customized near power boost**, tailored to the wearer's age and vision needs. Designed to **alleviate eyestrain** caused by frequent use of digital devices, these lenses provide optimal support when viewing digital screens.

0.50 D

Ideal wearer: those with mild symptoms of visual fatigue.


0.75 D

Ideal wearer: people with moderate symptoms of visual fatigue.

1.00 D

Ideal wearer: pre-presbyopes and those with severe symptoms of visual fatigue.

Excellent distance and
peripheral vision.



Concentrated power boost.



Fitting
Configuration

Fitting Configuration.

“ Fitted perfectly to you ”

→ Individual Fit

Optimizes the vision zones based on the unique data of the wearer, including their prescription, facial morphology, frame shape, and use position. If individual parameters are not supplied, default values will be used.

→ Frame Fit

Regardless of the shape or size of the frame, our technology adjusts the lens design structure to provide optimal performance.

→ Lens Fit

Provides a sleek and modern aesthetic for the wearer, even with a high prescription. This is done through utilizing the optimal base curves and thickness optimization.

Contact us

www.artisanlabnetwork.com/

Distributed
by:

