

Varilux®

Now available with
NEW Transitions® GEN S™

Varilux®

#1 progressive lens brand

recommended by optometrists and opticians¹



VARILUX® LENSES AVAILABILITY for Eyecare Professionals

Effective APRIL 2024

(1) Survey conducted in 2021 by an independent third party, sponsored by Essilor. Results were reported by independent opticians and optometrists.

(2) Crizal® coating is Essilor® lenses treatment that acts as an invisible shield of protection.



essilor

Artisan
LAB NETWORK

VARILUX® XR DESIGN

LENS ID



Product requires additional processing time.

Crizal

AR restricted to Crizal coatings.

MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	★ Clear	1.50									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.53	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN & LITE®	★ Clear	1.67									-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓						-12.00 to +9.00	0.75 to 4.00D	-6.00
	★ Clear**	1.74									-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™***	1.74	✓	✓							-14.00 to +13.00	0.75 to 4.00D	-6.00

VARILUX® XR TRACK

LENS ID



Product requires additional processing time.

Crizal

AR restricted to Crizal coatings.

MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	★ Clear	1.50									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.53	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN & LITE®	★ Clear	1.67									-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓						-12.00 to +9.00	0.75 to 4.00D	-6.00
	★ Clear**	1.74									-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™***	1.74	✓	✓							-14.00 to +13.00	0.75 to 4.00D	-6.00

MFH=14mm

NOTES:

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with Crizal.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ Currently Available

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



Artisan AR Treatments available.

Crizal

Crizal AR only available.

VARILUX® X DESIGN

LENS ID



Varilux X series all require additional processing time until September.

Crizal®

AR restricted to Crizal coatings until September.

MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	★ Clear	1.50									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.53	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN & LITE®	★ Clear	1.67									-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓						-12.00 to +9.00	0.75 to 4.00D	-6.00
	★ Clear**	1.74									-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™**	1.74	✓	✓							-14.00 to +13.00	0.75 to 4.00D	-6.00

VARILUX® X FIT

LENS ID



VARILUX® X 4D

LENS ID



Crizal®

AR restricted to Crizal coatings until September.

MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	ADD POWER*	MAX CYLINDER
PLASTIC	★ Clear	1.50									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.53	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +6.00	0.75 to 4.00D	-6.00
	Essilor Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓						-10.00 to +6.00	0.75 to 4.00D	-6.00
THIN & LITE®	★ Clear	1.67									-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-12.00 to +9.00	0.75 to 4.00D	-6.00
	Transitions® XTRActive® new generation	1.67	✓SM	✓	✓						-12.00 to +9.00	0.75 to 4.00D	-6.00
	★ Clear**	1.74									-14.00 to +13.00	0.75 to 4.00D	-6.00
	Transitions® GEN S™**	1.74	✓	✓							-14.00 to +13.00	0.75 to 4.00D	-6.00

MFH=14mm

NOTES:

*ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Preveincia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

★ Available in Essential Blue Series® lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ Currently Available

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



Artisan AR Treatments available.



Crizal AR only available.

VARILUX® PHYSIO® W3+  												VARILUX® PHYSIO® W3+ FIT  												PHYSIO® W3+ EYECODE™  											
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	BASE CURVES	ADD POWER*	CYLINDER																					
PLASTIC	★ Clear	1.50									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Essilor Polarized	1.50	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™	1.50	✓	✓	✓						-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™	1.53	✓	✓							-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® XTRActive® new generation	1.53	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25	0.75 to 3.50, 4.00D	-6.00																					
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Essilor Polarized	1.59	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +7.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25	0.75 to 3.50, 4.00D	-6.00																					
THIN & LITE®	★ Clear	1.60									-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™	1.60	✓	✓	✓						-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	★ Clear	1.67									-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Essilor Polarized	1.67	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Xperio UV™ Polarized	1.67	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓						-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00	0.75 to 3.50, 4.00D	-6.00																					
	Clear**	1.74									-18.00 to +13.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00																					
	Transitions® GEN S™**	1.74	✓	✓							-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00	0.75 to 3.50, 4.00D	-6.00																					

MFH=14mm

NOTE: *ADD Power in steps of 0.25D. At fitting height less than 17mm, ADD power is up to +3.50.

VARILUX® COMFORT MAX												VARILUX® COMFORT MAX FIT											
MATERIAL		PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE		BASE CURVES		ADD POWER*	CYLINDER						
PLASTIC	★ Clear	1.50										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Essilor Polarized	1.50	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.50	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Transitions® GEN S™	1.50	✓	✓	✓							-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Transitions® XTRActive® new generation	1.50	✓	✓	✓							-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
TREXA® (TRIVEX®)	★ Clear	1.53										-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Transitions® GEN S™	1.53	✓	✓	✓							-10.00 to +7.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
	Transitions® XTRActive® new generation	1.53	✓	✓								-10.00 to +6.00	2.00, 3.25, 4.50, 5.75, 7.25		0.75 to 3.50D	-6.00							
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59										-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D	-6.00							
	Essilor Polarized	1.59	✓	✓								-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.59	✓	✓								-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D	-6.00							
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D	-6.00							
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓							-10.00 to +6.00	1.75, 3.00, 4.00, 5.50, 7.25		0.75 to 3.50D	-6.00							
THIN & LITE®	★ Clear	1.60										-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00		0.75 to 3.50D	-6.00							
	Transitions® GEN S™	1.60	✓	✓	✓							-12.00 to +8.00	1.75, 3.00, 4.00, 5.50, 7.00, 8.00		0.75 to 3.50D	-6.00							
	★ Clear	1.67										-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D	-6.00							
	Essilor Polarized	1.67	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D	-6.00							
	Xperio UV™ Polarized	1.67	✓	✓								-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D	-6.00							
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.50D	-6.00							
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.50, 8.00		0.75 to 3.00D	-6.00							

VARILUX PHYSIO DRx™														LENS ID 		VARILUX COMFORT DRx™														LENS ID 	
MATERIAL	PRODUCT	INDEX	G	B	GG	SA	AMT	AMB	EM	RB	POWER RANGE	BASE CURVES					ADD POWER*		CYLINDER												
PLASTIC	★ Clear	1.50									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75					0.75 to 3.50, 4.00D [†]		-6.00												
	Essilor Polarized	1.50	✓	✓							-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25					0.75 to 3.50, 4.00D [†]		-6.00												
	Xperio UV™ Polarized	1.50	✓	✓							-10.00 to +7.50	2.00, 3.25, 4.25, 6.25, 7.25, 8.25					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® GEN S™	1.50	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® new generation	1.50	✓	✓	✓						-10.00 to +7.50	2.25, 4.25, 5.25, 8.25					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® Polarized™	1.50	✓									-10.00 to +7.50	2.00, 3.50, 4.25, 5.25, 6.25, 7.25, 8.75, 10.50					0.75 to 3.50, 4.00D [†]		-6.00											
TREXA® (TRIVEX®)	★ Clear	1.53									-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® GEN S™	1.53	✓	✓	✓						-10.00 to +6.00	2.00, 3.50, 4.25, 6.25, 7.25					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® new generation	1.53	✓	✓							-10.00 to +7.50	2.00, 3.50, 4.25, 6.25, 7.25, 8.75					0.75 to 3.50, 4.00D [†]		-6.00												
AIRWEAR® (POLYCARBONATE)	★ Clear	1.59									-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	Essilor Polarized	1.59	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	Xperio UV™ Polarized	1.59	✓	✓							-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® GEN S™	1.59	✓	✓	✓	✓	✓	✓	✓	✓	-10.00 to +6.00	1.50, 3.25, 4.25, 5.50, 7.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® new generation	1.59	✓SM	✓SM	✓						-10.00 to +8.00	1.50, 3.25, 4.25, 5.50, 7.00, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® Polarized™	1.59	✓									-10.00 to +8.00	1.50, 3.25, 4.25, 4.75, 5.50, 7.00, 8.50					0.75 to 3.50, 4.00D [†]		-6.00											
THIN & LITE®	★ Clear	1.60									-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® GEN S™	1.60	✓	✓							-12.00 to +8.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.50					0.75 to 3.50, 4.00D [†]		-6.00												
	★ Clear	1.67									-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Essilor Polarized	1.67	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Xperio UV™ Polarized	1.67	✓	✓							-14.00 to +9.00	1.75, 2.75, 3.75, 5.25, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® GEN S™	1.67	✓	✓	✓	✓	✓	✓	✓	✓	-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® new generation	1.67	✓SM	✓SM	✓						-14.00 to +9.00	1.75, 2.75, 3.75, 6.00, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00												
	Transitions® XTRActive® Polarized™	1.67	✓									-14.00 to +9.00	1.75, 2.75, 3.75, 4.75, 6.00, 6.75, 8.00					0.75 to 3.50, 4.00D [†]		-6.00											
	★ Clear**	1.74										-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00					0.75 to 3.50, 4.00D [†]		-6.00											
	Transitions® GEN S™**	1.74	✓	✓								-20.00 to +13.00	0.75, 1.75, 2.75, 3.75, 5.25, 6.50, 8.00, 10.00					0.75 to 3.50, 4.00D [†]		-6.00											

MFH=17mm

NOTES:

*ADD Power in steps of 0.25D.

**THIN & LITE 1.74 lenses come systematic with choice of Crizal® Prevencia®, Crizal Sapphire® 360° UV, Crizal Avancé UV™, Crizal® Rock™ or Crizal Alizé UV™.

***For the complete list of Transitions® color availability please refer to the *Transitions* Availability Chart.

For complete Sunwear availability please refer to the Essilor Sun Portfolio Availability Chart.

† 4.00D add only available in *Varilux Physio DRx*.

★ Available in *Essential Blue Series®* lenses.

G: Gray, B: Brown, GG: Graphite Green, SA: Sapphire, AMT: Amethyst, AMB: Amber, EM: Emerald, RB: Ruby

✓ **Currently Available**

SM: Available in Transitions® XTRActive® Style Mirrors in the following six flash mirror colors: Gold, Silver Shadow, Blue, Green, Pink, and Red.



Artisan AR Treatments available.



Crizal AR only available.

VARILUX PHYSIO SHORT DRx™

LENS ID

H 20 S

VARILUX COMFORT SHORT DRx™

LENS ID

H 20 S

RX AVAILABILITY FOR VARILUX PHYSIO SHORT DRX AND VARILUX COMFORT SHORT DRX MATCHES VARILUX PHYSIO DRX AND VARILUX COMFORT DRX EXCEPT ADD POWER RANGE IS **0.75 TO 3.00D** ACROSS ALL MATERIALS.

MFH=14mm

Varilux®

PROGRESSIVE LENSES

Sharp vision and smooth transitions at any distance

PERSONALIZED MEASUREMENTS

Required: ● Optional, Default Measurements Accepted: ■

	PD and Fitting Height	Pantoscopic Tilt	Wrap Angle	Vertex Distance	Near Vision Behavior	Eye Rotation Center Distance	Natural Head Position	Leading Dominant Eye	H/E Ratio	Stability Coefficient
Varilux® XR track	●	●	●	●	●					
Varilux® XR design	●	■	■	■						
Varilux® X 4D with H/E†	●	●	●			●	●	●	●	●
Varilux® X 4D†	●	●	●			●	●	●		
Varilux® X Fit	●	■	■	■						
Varilux® Physio® W3+ Eyecode™†	●	●	●			●	●			
Varilux® Physio® W3+ Fit	●	■	■	■						
Varilux® Comfort Max Fit	●	■	■	■						

† Varilux X 4D lenses, Varilux Physio W3+ Eyecode lenses are exclusive to the Visioffice® System.



(1) Survey conducted in 2021 by an independent third party, sponsored by Essilor. Results were reported by independent opticians and optometrists.

©2024 Essilor of America, Inc. All rights reserved. Unless indicated otherwise, all registered trademarks and trademarks are the property of Essilor International and/or its subsidiaries in the United States and in other countries. Transitions, Transitions Signature, XTRActive, and the Transitions logo are registered trademarks and is a trademark of Transitions Optical, Inc., used under license by Transitions Optical Ltd. GEN S is a trademark of Transitions Optical Limited. Photochromic performance is influenced by temperature, UV exposure, and lens material. These products may be protected by one or more patents listed at www.essilorusa.com/patents. 347550_PRO_VAR_SHK/HB 05/24



essilor.com